



Remote Control RC software application

The Remote Control (RC) software opens up the possibility for the user to fully control the Finapres® NOVA remotely. The user can keep distance from the subject while using a computer to control and monitor the measurement!

With the RC software, it is possible to set-up a connection with multiple devices:

- ✓ Multiple Finapres® NOVA devices can be linked to one computer
- ✓ Multiple computers can be linked to one Finapres® NOVA device

Requirements to use the Remote Control software

- ✓ Computer (operating on Windows) or tablet with the TightVNC software installed
- ✓ Ethernet network connection cable

Some examples of practical applications:

- ✓ X-ray research
- ✓ High-pressure chamber research
- ✓ Performing measurements on contagious patients (like COVID-19)
- ✓ Sleep (disorder) research
- ✓ Educational courses



Scan the QR-code to read more about the Remote Control software application on our website.