




NOVAscope PC software application

The Finapres® NOVA facilitates accurate review of recorded measurements!

Measurements files recorded with the Finapres® NOVA can be reviewed in two ways:

1. On the NOVA device itself, in Review Mode.
2. With the NOVAscope PC application*, installed on a separate computer or laptop.

* This application requires an additional license.



Measurement review with the Finapres® NOVA

Measurement review features

- ✓ Full measurement review
- ✓ Adding markers
- ✓ Adding measurement comments and print comments
- ✓ Create measurement reports for (digital) printing
- ✓ Beat-to-beat data export
- ✓ Raw data export
- ✓ GAT video synchronisation*

*This feature is only available with the NOVAscope PC application



Scan the QR-code to read more about the NOVAscope PC software application, on our website.